

TISHK INTERNATIONAL UNIVERSITY FACULTY OF APPLIED SCIENCE MEDICAL ANALYSIS DEPARTMENT



Autumn Semester 2023-2024 Course Name : **Biochemistry (Practical)** Stage : 3 Lecture 4: **Iodine Test** Lecture: Dr. Soma Majedi / Ph.D. in Organic Chemistry









Main Goal:

➢ Using Iodine to test for the presence of Starch is a common experiment.

Minor Goals:

1- Understanding somewhat about the structure of **Starch**.

2- Learning some more about the structure of Iodine-Starch Complex.

Principle:



A solution of elemental iodine (I_2) and potassium iodide (KI) in water (contain triiodide anion; I^{3-}) has a light **orange-brown** color (**Lugol's iodine**).

 $I_2 + I^- \rightarrow I_3^-$

If it is added to a sample that contains starch, produces an intense "**blue-black**" color. Starch forms a very **dark blue-black** complex with **triiodide**.



Such as the bread pictured above, the color changes to a **deep blue**



- 1. The complex is **not formed** if only iodine or only iodide (**I**-) is present.
- 2. The color of the starch complex is so deep, that it can be detected visually when the concentration of the iodine is as low as 0.00002 M at 20 °C.
- 3. The intensity of the color **decreases** with **increasing temperature** and with the presence of **water-miscible organic solvents** such as **ethanol**.





4. The test cannot be performed at very low pH (acidic) due to the Hydrolysis of the starch under these conditions.

5. Iodine by itself is very poorly soluble in water. One way to dissolve Iodine in water is to add **Potassium Iodide** or **Sodium Iodide**.

6. Foods that are high in starch include grain foods and some vegetables, such as dried beans and peas, potatoes, corn, and lentils.

7. As fruits **ripen**, the amount of starch can **decrease**.



Lugol's iodine is available in various concentrations.

The most commonly used 15% solution consists of: 10% (wt/v) Potassium Iodide (KI)

5% (wt/v) Iodine (I_2) mixed in distilled water has a total iodine content of 126.5 mg/mL.



The 15% solution thus has a total iodine content of 6.32 mg per drop of 0.05 mL; the 2% solution has 0.84 mg total iodine content per drop.

Starch is a carbohydrate found in plants. Natural starches - from plants consists of two different types of polysaccharides that are made up of glucose units which are connected in two different ways. One is the **linear amylose (10 - 20%)** and the other is the **branched amylopectin (80-90%)**. The **structure of the helical amylose** is key to the loding starsh reaction.

the Iodine-starch reaction. A helix is a coil or a spring.



Procedure of Iodine Test

Step1:

- 1- Pour **2 ml** of each sample solutions in a test tube.
- 2- Add 3 drops of Lugol's iodine (triiodide anion; I³⁻) to the solution.
 3- (Mix well).

Step 2:

- Record the observed color changes.

Note:

These conditions are whether sample used in the test is solid or a liquid.







Starch-Iodine-Complex





Dark blue-black, if the test is positive





+ result - result



In the case of presence of **starch**, the color of the solution changes.

The color change should be a **blueblack color** which appears only in the case if **positive results** which are that the **starch** is present in the solution.





Semi-Quantitative Test

Through visual comparisons or more advanced techniques like spectrophotometry to estimate approximate starch concentration.

Iodine Test for Starch





Interpretation:

Amylose is the compound that is responsible for the blue color. Its chain forms a helix shape, and iodine can be bound inside this helix.





Homework



- 1. Describe the "charge-transfer" and explain how this phenomenon affects on the iodine test?
- 2. Why does the blue color disappear after heating, or does it decrease in intensity by adding Ethanol?
- 3. Why does the amount of starch decrease as fruits ripen?
- ✓ <u>https://www.google.com/search?q=iodine+test+cartoon&source=lmns&tbm</u> =vid&bih=739&biw=1536&hl=en&sa=X&ved=2ahUKEwjqK08xJGCAxW2mycCHU85CVkQ_AUoAnoECAEQAg#fpstate=ive&vld =cid:1ee229c6,vid:6rjakS5W0Go,st:0
- ✓ <u>https://www.google.com/search?q=iodine+test+cartoon&source=lmns&tbm</u> =vid&bih=739&biw=1536&hl=en&sa=X&ved=2ahUKEwj-<u>qKO8xJGCAxW2mycCHU85CVkQ_AUoAnoECAEQAg#fpstate=ive&vld</u> =cid:3607c1a7,vid:HO_q8GPl3bE,st:0







References

Í.VA

?

Laboratory Manual for Practical Biochemistry, MK Ganesh, YM Shivaraja Shankara, AR Shivashankara, ISBN: 9789350902769, DOI: 10.5005/jp/books/12009_5, Edition 2/e, Publishing Year: 2013.

Practical Clinical Biochemistry: Methods and Interpretations, Chawla Ranjna, ISBN: 9789389188769, DOI: 10.5005/jp/books/18754, Edition: 5/e, Publishing Year: 2020.