



DETERMINATION OF SERUM UREA

Dr. Jaafaru Sani Mohammed/ Mr. Tolaz Khadhm

Advanced Clinical Biochemistry I

Fall Semester

Week Three

15/10/2024



Objectives

- At the end of this lab session, **students should understand:**
- The importance of Urea test in the diagnosis of kidney disease
- The principles behind the urea test
- The step-by-step procedures for urea test
- The clinical significance of urea tests


Introduction

- Ammonia can be produced by the breakdown of amino acids, or by the gut bacteria in humans.
- If the level of ammonia in the blood becomes too high, then it becomes toxic to the brain.
- The urea cycle removes ammonia from the blood and makes urea, which is eventually excreted as urine.
- This cycle is carried out by the cells of the liver.

Urea Cycle

- **Carbamoyl phosphate** is formed from NH_3 and bicarbonate, by carbamoyl phosphate synthetase (CPS) in the mitochondria.
- Ornithine transcarbamoylase (OTC) condenses carbamoyl phosphate and ornithine to form **citrulline** and then transported to the cytosol
- Argininosuccinate synthetase (AS) condenses citrulline and aspartate to form **argininosuccinate in the cytosol**.
- Argininosuccinate is broken down into **arginine** and fumarate by argininosuccinate lyase (AL).
- Arginine is broken down into **urea** and ornithine **by arginase**.
- **Ornithine translocase** transports ornithine into the mitochondria.

UREA CYCLE



Clinical Significance

- Elevated levels of urea can appear in blood (Uremia) in diets with excess protein.
- Renal disease, heart disease, gastrointestinal hemorrhage, dehydration, or renal obstruction could be developed due to defect in the cycle
- Clinical diagnosis should not be made on a single test result.

Principle

- Urea in the sample is hydrolyzed enzymatically into ammonium (NH_4^+) and CO_2 .
- Ammonium ions react with salicylate and hypochlorite (NaClO), in the presence of the catalyst nitroprusside, to form a blue-green indophenol as follows:



Reagents

- **R₁:** **Salicylate**
- **R₂:** **(Enzymes)** **Urease**
- **R₃:** **Alkaline reagent**
- **R₄:** **Standard (40mg/dL or 6.66 mmol/L)**

Procedure

Pipette into the tubes	Blank	Standard	Assay/Test
$R_W = R_1 + R_2$ (μL)	1000	1000	1000
dH ₂ O (μL)	5		
Standard (R_4) (μL)		5	
Serum (μL)			5
Mix the reagents and wait for 4 minutes at room temperature (T_R)			
R3 (Alkaline) (μL)	1000	1000	1000
Mix and allow to stand for 5 minutes at T_R . Read the absorbance at 600 nm against blank. The color of the reaction stays only for 2 minutes hence OD should be taken immediately.			

Calculation of result

○ Urea Conc. (mg/dl) = $\frac{A_{\text{sample}}}{A_{\text{standard}}} \times 40 \text{ (std conc.)}$

○ Conversion factor: $\text{mg/dl} \times 0.1665 = \text{mmol/L}$

○ Reference Value:

- Serum: 15 – 45 mg/dL (2.49 – 7.49 mmol/L)

Second Practical Report

- Refer to practical note 1 and write a step-by-step practical report for today's experiment.
- The report should comprise all the components outlined in the standard report format and should not be more than five (5) or less than three (3) pages.
- The submission date is Tuesday 22nd October 2024 during practical session

