



# **Demonstration of Energy Production During Cellular Respiration Using Yeast**

Course instructor: Jibril H. Yusuf PhD.

E-mail: [jibril.habib@tiu.edu.iq](mailto:jibril.habib@tiu.edu.iq)

Assistant: Ms. Lava Ali

Course: General Biology

Date 10-02-2026

# Outline

- Concept of Energy and Metabolism
- Cellular Respiration
- Why Yeast is Used in This Experiment.
- Principle of the Experiment
- Materials
- Procedure
- Expected Results

# ❖ **Concept of Energy and Metabolism:**




- **Energy:** The capacity to do work that bring about change living organism.
- **Metabolism:** The sum of all chemical reactions in a cell or organism.
- ✓ **Catabolism:** Breakdown of molecules to release energy
- ✓ **Anabolism:** Synthesis of molecules using energy

# ❖ Concept of Energy and Metabolism:




Feature	Catabolism	Anabolism
Meaning	Breaking down	Building up
Energy	Releases	Consumes
Molecules	Big → Small	Small → Big
Example	Respiration	Protein synthesis



# ❖ Concept of Energy and Metabolism:

- Cellular respiration is a catabolic process that releases energy stored in glucose.





# ❖ Cellular Respiration:



- Cellular respiration is the process by which cells convert glucose into ATP.

## ➤ General equation:



## ➤ In yeast, respiration can occur:

- Aerobically (with oxygen)
- Anaerobically (fermentation)



# ❖ Why Yeast is Used in This Experiment:



1. Yeast is a living organism with active metabolism.
2. Easy to handle and safe for teaching laboratories.
3. Rapidly metabolizes glucose.
4. Produces carbon dioxide, which can be easily observed

# ❖ Laboratory Experiment: Energy and Metabolism



➤ Demonstration of Energy Production During Cellular Respiration Using Yeast.

- **Aim:** To demonstrate energy production during metabolism by observing carbon dioxide ( $\text{CO}_2$ ) release during cellular respiration in yeast.

➤ **Principle:** All living cells require energy for metabolic activities. This energy is produced mainly in the form of ATP through cellular respiration.



➤ Yeast cells metabolize glucose using enzymatic reactions.

✓ During respiration:

- Glucose is broken down
- Energy is released
- Carbon dioxide (CO<sub>2</sub>) is produced as a by-product

➤ The production of CO<sub>2</sub> indicates active metabolism and energy release.

## ❖ **Materials Required:**

- Dry yeast
- Glucose solution (5–10%)
- Warm water (35–40°C)
- Test tubes
- Balloon or delivery tube
- Measuring cylinder
- Rubber band
- Test tube rack
- Stopwatch




## ❖ Procedure:




1. Label two test tubes as **Test** and **Control**.
2. Add **10 ml of warm water** to each test tube.
3. Add **1 teaspoon of yeast** to both tubes.
4. Add **5 ml glucose solution** to the **Test** tube only.
5. Do not add glucose to the **Control** tube.
6. Fix a balloon tightly over the mouth of each test tube using a rubber band.
7. Place the tubes in a warm area.
8. Observe the setup for **20–30 minutes**.

# ❖ Expected Result



## Energy Production During Metabolism: CO<sub>2</sub> Release in Yeast Respiration



- The balloon on the Test tube inflates due to CO<sub>2</sub> production.
- The balloon on the Control tube shows little or no inflation.

# ❖ Significance of the Experiment:



- Demonstrates that living cells produce energy.
- Shows the relationship between metabolism and energy.
- Helps understand basic biochemical processes relevant to:
  - Human physiology.
  - Medical laboratory analysis.



# References



- Reece, J. B., et al. (2014). Campbell Biology, 10th Ed. Pearson.
- Alberts, B. et al. (2015). Molecular Biology of the Cell, 6th Ed. Garland Science.
- Svitil, K. (2006). "Crenation and Osmosis in Cheek Cells". The American Biology Teacher, 68(6), 364–367.
- Oparka, K. J. (1994). Plasmolysis: New insights into an old process. New Phytologist, 126(4), 571–591. <https://doi.org/10.1111/j.1469-8137.1994.tb02958.x>





**Thanks**