

Tishk International University
Faculty of Applied Science
Nutrition and Dietetics Department
Spring semester
2nd Grade



Nutritional Biochemistry II

Qualitative Determination of Enzymes

Lecturer: Amani Tahsin





Objectives

- Qualitative enzyme tests provide a simple and effective method for detecting enzyme activity through observable physical or chemical changes.
- These methods are widely used in clinical, microbiological, and biochemical laboratories for preliminary enzyme identification.



Qualitative determination of enzymes

- Refers to identifying the *presence* and *activity* of a specific enzyme in a sample without measuring its exact concentration or activity rate quantitatively.
- In biochemistry practical courses, qualitative tests are commonly based on:
 1. **Observable color change**
 2. **Precipitate formation**
 3. **Gas production**
 4. **Substrate disappearance**
 5. **Product formation**



General Principle

Enzymes catalyze specific biochemical reactions.

If the enzyme is present and active → it converts a substrate into a product → a visible change occurs.




Common Qualitative Enzyme Tests

1. Catalase Test

Principle:

Catalase breaks down hydrogen peroxide into water and oxygen.



Procedure:

Add a few drops of hydrogen peroxide (H_2O_2) to the sample.

Observe immediately.

Positive result:

✓ Rapid bubble formation (oxygen gas)

Negative result:

✗ No bubbling



2. Amylase Test (Iodine Test)

Principle:

Amylase hydrolyzes starch into maltose.

Iodine gives a blue-black color with starch.

Procedure:

Incubate starch with enzyme sample.

Add iodine solution.

Positive result (amylase present):

✓ No blue-black color (starch digested)

Negative result:

✗ Blue-black color remains



3. Urease Test

Principle:

Urease breaks down urea into ammonia and CO_2 .
Ammonia increases pH.

Procedure:

Add sample to urea broth containing phenol red.

Positive result:

✓ Pink/red color (alkaline)

Negative result:

✗ Yellow/orange (no pH change)



4. Protease Test

Principle:

Protease digests proteins (e.g., casein in milk agar).

Procedure:

Inoculate sample on milk agar.

Incubate.

Positive result:

✓ Clear zone around colony (protein hydrolysis)

Negative result:

✗ No clearing



5. Lipase Test

Principle:

Lipase hydrolyzes fats into fatty acids and glycerol.

Observation methods:

Formation of clear zones on tributyrin agar

Color change using pH indicators



Important Factors Affecting Qualitative Tests

- Temperature
- pH
- Incubation time
- Enzyme concentration
- Substrate specificity



Precautions

- Maintain optimal temperature (37°C unless stated otherwise).
- Use fresh reagents.
- Avoid contamination of samples.
- Handle hydrogen peroxide carefully.
- Record observations immediately.