



Candida Species Lab Diagnosis

Course: Medical Mycology
Week 2

Lab Date 27/8/ 2026



Mr. Soran Kayfi & Ms. Bana Wali

Pathogenic Scope

Candida species are opportunistic fungal pathogens residing as normal microbiota on human mucosal surfaces.

Under immunocompromised conditions, they cause **thrush**, **vaginitis**, and **invasive candidiasis infections**.

Primary Pathogen: *Candida albicans* accounts for 50–70% of infections, though non-*albicans* species are emerging rapidly.



thrush

Diagnosis of *Candida* spp. Infections

- The diagnosis of *Candida* infections depends on the type of infection (e.g., superficial vs. invasive) and the clinical setting.
- *Candida albicans*: The most prevalent and well-known species, often responsible for infections.

The common diagnostic method

- Clinical Examination.
- Laboratory Tests.
- Antifungal Susceptibility Testing.

Clinical Examination

Oral Candidiasis: White patches on the tongue, inner cheeks, and throat.

Cutaneous Candidiasis: Red, itchy rash in warm, moist areas (e.g., groin, under breasts).

Invasive Candidiasis: *Fever, chills, and organ dysfunction, particularly in immunocompromised patients.*



Oral Candidiasis



Skin Candidiasis



Invasive Candidiasis

2. Laboratory Tests

- A. Microscopy & Staining.
- B. Culture Methods.
- C. Molecular & Serological Tests.
- D. Histopathology (for Tissue Infections).

A. Microscopy & Staining

- **Potassium hydroxide (KOH) wet mount:** Detects yeast cells and pseudohyphae in skin, vaginal, or oral samples. **Gram stain:** Shows Gram-positive budding yeast.
- **Calcofluor white stain: Fluorescent stain** that binds to fungal cell walls for better visualization.

Microscopy & Staining

- **KOH wet mount:** Detects **yeast cells** and **pseudohyphae** in skin, vaginal, or oral samples, it is a diagnostic tool used to identify fungal infections, particularly those caused by *Candida albicans*.

Sample Collection

- **Oral candidiasis:** Swab from the tongue, buccal mucosa, or throat.
- **Vaginal candidiasis:** Vaginal swab or discharge sample.
- **Cutaneous candidiasis:** Skin scrapings from affected areas.
- **Nail candidiasis:** Nail scrapings or clippings.

Main steps of KOH Wet Mount



- **Preparation:** A small amount of skin scrapings, vaginal discharge, or oral swab is collected from the affected area. The sample is placed on a glass slide.
- **Addition of KOH:** A drop of 10%–20% KOH is added to the sample, let the slide sit for 5–10 minutes at room temperature to dissolve keratin and cellular debris, making fungal elements more visible.
- **Examination:** The mixture is covered with a coverslip and examined under a microscope.

Microscopic Examination

- Use a light microscope under low (10x) and high (40x) magnification.
- Look for budding yeast cells and pseudohyphae, which are characteristic of *Candida albicans*.




Interpretation of Results

Finding

Budding yeast cells

Pseudohyphae (elongated yeast cells)

No fungal elements




Interpretation

Suggests *Candida* species, including *C. albicans*.


Highly suggestive of *Candida albicans*.

Negative result, but does not rule out infection (culture may be needed).




Budding yeast cells


Microscopic Examination



Hyphae and spores of *Candida albicans*





Microscopic Examination of *C. albicans* using gram stain



B. Culture Methods

1. **SDA:** Supports fungal growth; colonies appear creamy-white.
2. **CHROMagar Candida:** Differentiates Candida species by colony color.
3. **Blood Culture:** Detects Candida in bloodstream infections.



SDA



CHROMagar

Candida albicans on SDA

- SDA is a selective medium that supports the growth of fungi, including yeasts like *Candida albicans*.

Colony Morphology

When cultured on SDA, *Candida albicans* typically exhibits the following characteristics:

- **Colony Color:** Colonies are usually white to cream-colored.
- **Colony Texture:** The texture can be smooth or slightly wrinkled, with a moist appearance.
- **Colony Size:** Colonies can grow to a diameter of 2-5 mm within 24-48 hours at 30-37°C.

Growth Conditions

- **Temperature:** Optimal growth occurs at 30-37°C.
- **pH:** SDA has a pH range of 5.6, which is conducive for fungal growth.
- **Oxygen Requirement:** *Candida albicans* is a **facultative anaerobe**, meaning it can grow in **both aerobic** and **anaerobic** conditions.

***Candida albicans* on Chromagar**

Chromagar is particularly useful in clinical microbiology for the rapid identification of *Candida* species, allowing for timely diagnosis and treatment of infections caused by these fungi.

Candida albicans on Chromagar

- Chromagar Composition: Chromagar is a **selective medium** that contains various nutrients and chromogenic substrates that allow for the differentiation of yeast species based on their metabolic activity.

The key **components** include:

1. **Peptones:** Provide nitrogen and essential growth factors. Chromogenic
2. **Substrates:** Allow for color development specific to certain *Candida* species.
3. **Agar:** Serves as a solidifying agent.

Candida albicans on Chromagar

- When *Candida albicans* is cultured on Chromagar, it typically produces a characteristic color change.
- A sample suspected of containing *Candida albicans* is streaked onto the Chromagar plate.
- **Incubation:** The plate is incubated at 35-37°C for 24-48 hours.
- **Observation:** After incubation, colonies of *Candida albicans* will appear as greenish-blue colonies due to the hydrolysis of the chromogenic substrate.


Candida albicans on Chromagar

Interpretation of Results


- **Greenish-blue colonies:** Indicative of *Candida albicans* presence.
- **Other colors:** Different *Candida* species may produce various colors, aiding in the differentiation of other species such as *Candida tropicalis* (blue colonies) and *Candida krusei* (pink colonies).

Candida albicans on Chromagar

Interpretation of Results



Candida albicans from different clinical samples



Culture and colony morphology of C. albicans.

